



Pharmacy Information System (PhIS) and Clinic Pharmacy System (CPS)

Klinik Kesihatan Bukit Kuda

**JKN SELANGOR
(ZONE 10)**

1ST Preventive Maintenance Report

P.M. Date	20 FEBRUARY 2018
Report ID	PM/REPORT/ZONE10/KKBUKITKUDA
Report Version	V1.0



PHIS & CPS PROJECT



© 2011 – 2018 Pharmaniaga Logistics Sdn Bhd (PLSB)

CONFIDENTIAL COPYRIGHTED MATERIAL – *The information includes all concepts, comments, recommendations, and material, contained herein shall remain the property of Pharmaniaga Logistics Sdn Bhd (PLSB). No portion of this document shall be disclosed, duplicated or used in whole or in part of any purpose other than the purpose of the Pharmacy Information System & Clinic Pharmacy System (PhIS & CPS) Project execution only.*



PHIS & CPS PROJECT



Revision History:

Ver. No	Date of Release	Consolidated by	List of Changes from Previous Version	Signature & Date
1.0	30/02/2018	Muhammad Anuar bin Awaluddin	Initial document	 MUHAMMAD ANUAR B. AWALUDDIN PHIS & CPS Implementer Pharmaliaga Logistik Sdn. Bhd. 22/2/18

Verified and Acknowledged by facility:

Received / Acknowledged by facility	
Name	
Signature / Company Stamp	 HANUM M. ABD RAHMAN Pegawai Farmasi U44 25/2/18
Designation / Department	
Date	



PhIS & CPS PROJECT



Appendix A

Certificate of Acceptance



EXECUTION OF PREVENTIVE MAINTENANCE

Project: PhIS & CPS

**1ST Preventive Maintenance (2018)
For
Klinik Kesihatan Bukit Kuda**

Summary findings for Preventive Maintenance conducted on 20/02/2018

No	Item	Allocated Quantity	Serviced Quantity	Remarks
1	Desktop	1	1	
2	LaserJet Printer	1	1	
3	Switch	1	1	
4	Server	1	1	
5	UPS	1	1	

Signature & Stamp:

<p><u>For PLSB</u></p>  Name: MUHAMMAD ANUAR D. AWALUDDIN Designation: PhIS & CPS Implementer Date: Pharmaniaga Logistic Sdn. Bhd.	<p><u>For Facility</u></p>  Name: HANUM WATI NUR ABD RAHMAN Designation: Pegawai Farmasi U44 Date: 20/2/18
--	--



PhIS & CPS PROJECT



Appendix B

Preventive Maintenance Checklist

**PREVENTIVE MAINTENANCE CHECKLIST
(SERVER / UPS / STORAGE)**

(FORM 2)

FACILITY : **KLINIK KESIHATAN BUKIT KUDA**

DEPT : **IT**

LOCATION : **Server Room**

STATE : **SELANGOR**

ZONE : **10**

SERVER 1	SERVER 2	UPS
PAT ID : CPS/SEL-03/KKS-03/SVR-001	PAT ID : CPS/SEL-03/KKS-03/UPS-001	
SERIAL NO : 06GH-CPC	SERIAL NO : 1517403382AN413	
MODEL : IBM X3500 M5	MODEL : LIEBERT GXT3-1000RT230	

! Alerts tick (V), fill in & elaborate details in FINDINGS / COMMENTS / DETECTION session.

PART I - HARDWARE , CABLING & CONNECTIVITY (INSPECTION & PHYSICAL CHECK)

Standard Task:

Cleaning

Verify Condition :

- Connection via LAN / USB
 - Fan/ Power Supply
 - Light Path Diagnosis panel
 - UTP Cable connect to switch
 - HA Gateway to Core Switch
- (For Hospital only)*

Function	SERVER 1	
	Not Function & Resolved	Net Function & Unresolved
<input checked="" type="checkbox"/>	<input type="checkbox"/>	Alert <input type="checkbox"/> Ticket No
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Function	SERVER 2	
	Not Function & Resolved	Net Function & Unresolved
<input type="checkbox"/>	<input type="checkbox"/>	Alert <input type="checkbox"/> Ticket No
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Function	UPS	
	Not Function & Resolved	Net Function & Unresolved
<input checked="" type="checkbox"/>	<input type="checkbox"/>	Alert <input type="checkbox"/> Ticket No
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

- By pass Test
- Webboard
- Multilink

Available	Type	Faulty
<input type="checkbox"/>		<input type="checkbox"/>

Qty

PART II - SETUP & HEALTH CHECK

Function	SERVER 1		SERVER 2	
	Not Function & Resolved	Alert	Not Function & Unresolved	Ticket No
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Run DSA & collect DSA Log

SERVER 1

SERVER 2

Software/Firmware Version :

Current	Updated
ClamAV 0.98.3	
Zabbix 2.2.6	
Bacula 5.2.12	
Glassfish 3.1.2.2	
UEFI	
iMM2	

Current	Updated

System Software

SUSE LINUX (OS)	
Kernel	
GCC	

Housekeeping:

- Clear /tmp folder > 2 weeks
- Archive Log Files

SERVER 1	
Current	Updated
SUSE 11	
3-0.76-0.11	
4.3.4	

SERVER 1	
Current	Updated

SERVER 2	
Current	Updated

SERVER 2	
Current	Updated

N/A

**PREVENTIVE MAINTENANCE CHECKLIST
(SERVER / UPS / STORAGE)**

(FORM 2)

PAT ID : _____
 SERIAL NO: _____
 TYPE : STORAGE

N/A

Function	Net Function & Resolved		Net Function & Unresolved	
	Alert	Ticks	Alert	Ticks
Perform Health Check:	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Storage Manager Error Alerts	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Storage events logs	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
HDD indicator lights	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Enable read/write cache	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Enable Media Scan both at controller level and LUN level	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Software/Firmware Version :

Controller	Updated
ESM	
Driver	

PART IV - FINDINGS / ALERTS / DETECTION

FURTHER ACTION / RECOMMENDATION

PART V - ACKNOWLEDGEMENT

CUSTOMER

Signature & Stamp
 Customer Name: **HANUM MAISYAH RAHMAN**
 Pegawai Masani U44
 Date: 20/07/18

PHARMANIAGA LOGISTICS SDN BHD (PLSB)

Signature & Stamp
 Engineer Name: **MUHAMMAD ANUAR B. AWALUDIN**
 Date: 20/03/18
 PHS & CPS Implementer
 Pharamiaga Logistic Sdn. Bhd.

**PREVENTIVE MAINTENANCE CHECKLIST
(CORE SWITCH / DISTRIBUTION SWITCH)**

(FORM 3)

FACILITY : **KLINIK KESIHATAN BUKIT KUDA** DEPT : **IT** STATE : **SELANGOR** ZONE : **10**

LOCATION :

CORE SWITCH / DISTRIBUTION SWITCH

PAT ID : **CPS/SEL-03/KKS-03/DSW-001**

SERIAL NO : **DN1190709VA**

Location : **Server Room**

Loaned Unlocated Relocated

PART I - STANDARD TASK

Cleaning

Verify Condition :
Cable & Connection
Power Supply
Check LED Status
System LED
Master LED
STACK ID

Function	Not Function & Resolved	Alert	Not Function & Unresolved	Ticket No
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

PART III - STANDARD SETUP

SSH Access & Connection
Execute Putty to start
Specify the IP to connect
Test Login
Test & Verify Management IP via console access

Function	Not Function & Resolved	Alert	Not Function & Unresolved	Ticket No
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

CORE SWITCH / DISTRIBUTION SWITCH

PAT ID : _____

SERIAL NO : _____

Location : _____

Loaned Unlocated Relocated

PART I - STANDARD TASK

Cleaning

Verify Condition :
Cable & Connection
Power Supply
Check LED Status
System LED
Master LED
STACK ID

Function	Not Function & Resolved	Alert	Not Function & Unresolved	Ticket No
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

PART III - STANDARD SETUP

SSH Access & Connection
Execute Putty to start
Specify the IP to connect
Test Login
Test & Verify Management IP via console access

Function	Not Function & Resolved	Alert	Not Function & Unresolved	Ticket No
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

PART IV - FINDINGS / COMMENTS / DETECTION

FURTHER ACTION / RECOMMENDATION

PART V - ACKNOWLEDGEMENT

CUSTOMER

Signature & Stamp
Customer Name
Date

HANU MASARAH ABD RAHMAN
Pharmaniaga Logistic Sdn Bhd
11/1/2022

PHARMANIAGA LOGISTICS SDN BHD (PLSB)
Signature & Stamp
Engineer Name
Date

MUHAMMAD AKUAP B. AWALLOON
POS & GPS Implementer
Pharmaniaga Logistic Sdn Bhd
28/1/22

PREVENTIVE MAINTENANCE CHECKLIST
(DEVICE)

FACILITY : **KLINIK KESIHATAN BUKIT KUDA** DEPT : **OUTPATIENT PHARMACY** LOCATION : **STATE : SELANGOR** ZONE : **10**

PAT ID : **CPS/SEL-08/KS-03/PCD-001** DEPT : **OUTPATIENT PHARMACY** LOCATION : **OUTPATIENT PHARMACY**

S/N CPU : **PC05FKXS** MONITOR : **V901R4GA** KEYBOARD : **0001208** MOUSE : **H5516H415N9** Leaned Unleased Relocated

PART I - STANDARD TASK		PART II - STANDARD HEALTH CHECK		PART III - STANDARD SETUP	
Function	Verify Condition :	Function	Alert	Function	Not Function & Resolved
All Device: Cleaning	<input checked="" type="checkbox"/>	Connection via LAN / USB	<input type="checkbox"/>	Auto update disabled:	<input type="checkbox"/>
Desktop/Notebook:	Yes <input type="checkbox"/> No <input checked="" type="checkbox"/> N/A <input type="checkbox"/>	Fan/ Power Supply/ Charger	<input type="checkbox"/>	(a) Java (JDK 6 update 30)	<input type="checkbox"/>
Device is harden	<input checked="" type="checkbox"/>	Hard Disk Utilization (GB) :	Free <input type="checkbox"/> Alert <input type="checkbox"/> Ticket No	(b) Internet Explorer 9 (64bits)	<input type="checkbox"/>
MEMORY :	2 GB <input checked="" type="checkbox"/> 4 GB <input type="checkbox"/>	1st Partition	<input type="checkbox"/>	CiamAv - Ver 0.98.3	<input type="checkbox"/>
		2nd Partition	<input type="checkbox"/>	Microsoft Office	<input type="checkbox"/>
			<input type="checkbox"/>	BCP (BCP PC only)	<input type="checkbox"/>
			<input type="checkbox"/>	Backup PC (Bacula Web)	<input type="checkbox"/>

PAT ID : **CPS/SEL-08/KS-08/LJP-001** DEPT : **OUTPATIENT PHARMACY** LOCATION : **OUTPATIENT PHARMACY**

S/N CPU : **451457LMZGM1Y** MONITOR : **MOUSE:** Leaned Unleased Relocated

PART I - STANDARD TASK		PART II - STANDARD HEALTH CHECK		PART III - STANDARD SETUP	
Function	Verify Condition :	Function	Alert	Function	Not Function & Resolved
All Device: Cleaning	<input checked="" type="checkbox"/>	Connection via LAN / USB	<input type="checkbox"/>	Auto update disabled:	<input type="checkbox"/>
Desktop/Notebook:	Yes <input type="checkbox"/> No <input checked="" type="checkbox"/> N/A <input type="checkbox"/>	Fan/ Power Supply/ Charger	<input type="checkbox"/>	(a) Java (JDK 6 update 30)	<input type="checkbox"/>
Device is harden	<input type="checkbox"/>	Hard Disk Utilization (GB) :	Free <input type="checkbox"/> Alert <input type="checkbox"/> Ticket No	(b) Internet Explorer 9 (64bits)	<input type="checkbox"/>
MEMORY :	2 GB <input type="checkbox"/> 4 GB <input type="checkbox"/>	1st Partition	<input type="checkbox"/>	CiamAv - Ver 0.98.3	<input type="checkbox"/>
		2nd Partition	<input type="checkbox"/>	Microsoft Office	<input type="checkbox"/>
			<input type="checkbox"/>	BCP (BCP PC only)	<input type="checkbox"/>
			<input type="checkbox"/>	Backup PC (Bacula Web)	<input type="checkbox"/>

PART IV - FINDINGS / ALERTS / DETECTION

FURTHER ACTION / RECOMMENDATION

PART V - ACKNOWLEDGEMENT	
CUSTOMER	PHARMANIAGA LOGISTICS SDN BHD (PLS)
Signature & Stamp	Signature & Stamp
Customer Name	Customer Name
Date	Date

(Handwritten Signature)
HANUZZA SARAH ABD RAHMAN
Pharmaniaga Logistik Sdn. Bhd.

(Handwritten Signature)
MUHAMMAD ANUAR B. AWALUDIN
PhIS & CPS Implementation
Pharmaniaga Logistic Sdn. Bhd.



PHIS & CPS PROJECT



Appendix C

Hardware Inventory List



PHIS & CPS PROJECT



Hardware Inventory List

List of equipment available in the facility as per below table:

1. PCD

No.	CPU	MONITOR	KEYBOARD	MOUSE	ASSET TAGGING	DEPARTMENT	LOCATION
1	PC05FKXS	V901R4GA	0001208	HS516HA15N9	CPS/SEL-03/KKS-03/PCD-001	Outpatient Pharmacy	Outpatient Pharmacy

2. LASERJET PRINTER

No.	LASERJET	ASSET TAGGING	DEPARTMENT	LOCATION
1	451457LM2GM1Y	CPS/SEL-03/KKS-03/LJP-001	Outpatient Pharmacy	Outpatient Pharmacy

3. SWITCH

No.	SERVER	ASSET TAGGING	DEPARTMENT	LOCATION
1	DNI190709VA	CPS/SEL-03/KKS-03/DSW-001	Outpatient Pharmacy	Server Room

4. SERVER

No.	SERVER	ASSET TAGGING	DEPARTMENT	LOCATION
1	06GHCPC	CPS/SEL-03/KKS-03/SVR-001	Outpatient Pharmacy	Server Room

5. UPS

No.	UPS	ASSET TAGGING	DEPARTMENT	LOCATION
1	1517403382AN413	CPS/SEL-03/KKS-03/LPS-001	Outpatient Pharmacy	Server Room