



Pharmacy Information System (PHIS) and Clinic Pharmacy System (CPS)

Klinik Kesihatan Kuala Perlis

**JKN PERLIS
(ZONE 1)**

1st Preventive Maintenance Report

P.M. Date	28 FEBRUARY 2019
Report ID	PM/REPORT/ZONE1/KKKUALAPERLIS
Report Version	V1.0



PHIS & CPS PROJECT

pharmaniaga

© 2011 – 2019 Pharmaniaga Logistics Sdn Bhd (PLSB)

CONFIDENTIAL COPYRIGHTED MATERIAL – *The information includes all concepts, comments, recommendations, and material, contained herein shall remain the property of Pharmaniaga Logistics Sdn Bhd (PLSB). No portion of this document shall be disclosed, duplicated or used in whole or in part of any purpose other than the purpose of the Pharmacy Information System & Clinic Pharmacy System (PhIS & CPS) Project execution only*

**Revision History:**

Ver. No	Date of Release	Consolidated by	List of Changes from Previous Version	Signature & Date
1.0	28/02/2019	Muhamad Redzuan Bin Aziz	Initial document	 Muhamad Redzuan Bin Aziz Field Technical Engineer 28/02/2019

Verified and Acknowledged by facility:

Received / Acknowledged by facility	
Name	
Signature / Company Stamp	 NOOR SHAHRUL BIN MUSA RPh 6411 Pegawai Farmasi Klinik Kesihatan Kuala Perlis
Designation / Department	
Date	28/2/19



PHIS & CPS PROJECT

pharmaniaga

Appendix A

Certificate of Acceptance

EXECUTION OF PREVENTIVE MAINTENANCE

Project: PhIS & CPS

**1st Preventive Maintenance (2019)
For
Klinik Kesihatan Kuala Perlis**

Summary findings for Preventive Maintenance conducted on 28/02/2019

No	Item	Allocated Quantity	Serviced Quantity	Remarks
1	Desktop	3	3	
2	LaserJet Printer	1	1	
3	Switch	1	1	
4	Server	1	1	
5	UPS	1	1	
6	Thermal Printer	1	1	

Signature & Stamp:

<p><u>For PLSB</u></p> <div style="text-align: right; margin-top: 20px;">  Muzhammad Redzuan Bin Aziz Fluid Technical Engineer </div> <p>Name: _____ Designation: _____ Date: 28/02/2019</p>	<p><u>For Facility</u></p> <div style="text-align: right; margin-top: 20px;">  NOOR SHAHRUL BIN MUSA RPh 6411 Pegawai Farmasi Klinik Kesihatan Kuala Perlis </div> <p>Name: _____ Designation: _____ Date: 28/2/19</p>
--	--



Appendix B

Preventive Maintenance Checklist

PREVENTIVE MAINTENANCE CHECKLIST (DEVICE)

FACILITY : **KLINIK KESIHATAN KUALA PERLIS** DEPT : **OUTPATIENT PHARMACY** LOCATION : **PERLIS** STATE : **PERLIS** ZONE : **1**

PAT ID : **CPS/PLS-01/KKS-05/PCD-001** DEPARTMENT : **DESKTOP / NOTEBOOK / LASERJET / THERMAL PRINTER** Location : **Screening, Verification, Labelling & Dispensing**

S/N CPU : **PC05FLR** MONITOR : **V901Z14** MOUSE : **H5515HA1K2Z** Keyboard : **0002536** Relocated : Unlocated : Loaned :

PART I - STANDARD TASK N/A

Verify Condition :
 Connection via LAN / USB
 Fan/ Power Supply/ Charger
 Hard Disk Utilization (GB):
 1st Partition
 2nd Partition

Desktop/Notebook: Yes No N/A
 Device is harden
 MEMORY : 2 GB
 4 GB

PART II - STANDARD HEALTH CHECK N/A

Function	Not Function & Resolved	Alert	Ticket No
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Used	Free	Alert	Ticket No
119	119	<input type="checkbox"/>	
76	136	<input type="checkbox"/>	

Verify Condition:
 IWP
 KYM
 PHIS Portal

Auto update disabled:
 (a) Java (JDK 6 update 30)
 (b) Internet Explorer 9 (64bits)
 ClamAv - Ver 0.99.1
 Microsoft Office
 BCP (BCP PC only)
 Backup PC (Bacula Web)

PAT ID : **CPS/PLS-01/KKS-05/PCD-002** DEPARTMENT : **DESKTOP / NOTEBOOK / LASERJET / THERMAL PRINTER** Location : **Screening, Verification, Labelling & Dispensing 1**

S/N CPU : **PC04MX7U** MONITOR : **V9022NMW** MOUSE : **H5509HA0K38** Keyboard : **01282950** Relocated : Unlocated : Loaned :

PART I - STANDARD TASK N/A

Verify Condition :
 Connection via LAN / USB
 Fan/ Power Supply/ Charger
 Hard Disk Utilization (GB):
 1st Partition
 2nd Partition

Desktop/Notebook: Yes No N/A
 Device is harden
 MEMORY : 2 GB
 4 GB

PART II - STANDARD HEALTH CHECK N/A

Function	Not Function & Resolved	Alert	Ticket No
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Used	Free	Alert	Ticket No
121	117	<input type="checkbox"/>	
76	136	<input type="checkbox"/>	

Verify Condition:
 IWP
 KYM
 PHIS Portal

Auto update disabled:
 (a) Java (JDK 6 update 30)
 (b) Internet Explorer 9 (64bits)
 ClamAv - Ver 0.99.1
 Microsoft Office
 BCP (BCP PC only)
 Backup PC (Bacula Web)

PAT ID : **CPS/PLS-01/KKS-05/PCD-003** DEPARTMENT : **DESKTOP / NOTEBOOK / LASERJET / THERMAL PRINTER** Location : **Medication Counselling & MTAC**

S/N CPU : **PC05FLR** MONITOR : **V902D1F8** MOUSE : **44F5738** Keyboard : **0005120** Relocated : Unlocated : Loaned :

PART I - STANDARD TASK N/A

Verify Condition :
 Connection via LAN / USB
 Fan/ Power Supply/ Charger
 Hard Disk Utilization (GB):
 1st Partition
 2nd Partition

Desktop/Notebook: Yes No N/A
 Device is harden
 MEMORY : 2 GB
 4 GB

PART II - STANDARD HEALTH CHECK N/A

Function	Not Function & Resolved	Alert	Ticket No
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Used	Free	Alert	Ticket No
121	117	<input type="checkbox"/>	
76	136	<input type="checkbox"/>	

Verify Condition:
 IWP
 KYM
 PHIS Portal

Auto update disabled:
 (a) Java (JDK 6 update 30)
 (b) Internet Explorer 9 (64bits)
 ClamAv - Ver 0.99.1
 Microsoft Office
 BCP (BCP PC only)
 Backup PC (Bacula Web)

Function	Not Function & Resolved	Current Version	Updated Version
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

**PREVENTIVE MAINTENANCE CHECKLIST
(DEVICE)**

FACILITY: **KLINIK KESIHATAN KUALA PERLIS** DEPT: **OUTPATIENT PHARMACY** LOCATION: **PERLIS** STATE: **PERLIS** ZONE: **1**

PAT ID : **CPS/PLS-01/KKS-05/LJP-001** DEPT: **OUTPATIENT PHARMACY** LOCATION: **PERLIS** STATE: **PERLIS** ZONE: **1**

S/N CPU : **45146PLM2X2BR** MONITOR: **N/A** KEYBOARD: **N/A** MOUSE: **DESKTOP / NOTEBOOK / LASERJET / THERMAL-PRINTER**

PART I - STANDARD TASK

All Device: Cleaning Yes No N/A

Desktop/Notebook: Device is harden Yes No N/A

MEMORY : 2 GB 4 GB

PART II - STANDARD HEALTH CHECK

Function	Not Function & Resolved	Alert	Ticket No
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Verify Condition :
 Connection via LAN / USB
 Fan/ Power Supply/ Charger
 Hard Disk Utilization (GB):
 1st Partition
 2nd Partition

Verification:
 IWP
 KYM
 PHIS Portal

Auto update disabled:
 (a) Java (JDK 6 update 30)
 (b) Internet Explorer 9 (64bits)
 ClamAv - Ver 0.99.1
 Microsoft Office
 BCP (BCP PC only)
 Backup PC (Bacula Web)

PAT ID : **CPS/PLS-01/KKS-05/THP-001** DEPT: **OUTPATIENT PHARMACY** LOCATION: **PERLIS** STATE: **PERLIS** ZONE: **1**

S/N CPU : **2305Y001491** MONITOR: **N/A** KEYBOARD: **N/A** MOUSE: **DESKTOP / NOTEBOOK / LASERJET / THERMAL PRINTER**

PART I - STANDARD TASK

All Device: Cleaning Yes No N/A

Desktop/Notebook: Device is harden Yes No N/A

MEMORY : 2 GB 4 GB

PART II - STANDARD HEALTH CHECK

Function	Not Function & Resolved	Alert	Ticket No
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Verify Condition :
 Connection via LAN / USB
 Fan/ Power Supply/ Charger
 Hard Disk Utilization (GB):
 1st Partition
 2nd Partition

Verification:
 IWP
 KYM
 PHIS Portal

Auto update disabled:
 (a) Java (JDK 6 update 30)
 (b) Internet Explorer 9 (64bits)
 ClamAv - Ver 0.99.1
 Microsoft Office
 BCP (BCP PC only)
 Backup PC (Bacula Web)

PART IV - FINDINGS / ALERTS / DETECTION

FURTHER ACTION / RECOMMENDATION

PART V - ACKNOWLEDGEMENT

CUSTOMER

Signature & Stamp: **NOOR SHAHRUL BIN MUSA**
 RPh 6411
 Pegawai Farmasi
 Klinik Kesihatan Kuala Perlis

Signature & Stamp: **Muhammad Redzuan Bin Aziz**
 Field Technical Engineer

Date: **28/12/19** Date: **28/12/2016**

**PREVENTIVE MAINTENANCE CHECKLIST
(SERVER / UPS / STORAGE)**

(FORM 2)

FACILITY : **KLINIK KESIHATAN KUALA PERLIS**

DEPT : **IT**

LOCATION : **Outpatient Pharmacy**

STATE : **PERLIS**

ZONE : **1**

SERVER 1	SERVER 2	UPS
PAT ID : CPS/PLS-01/KKS-05/SVR-001	PAT ID : CPS/PLS-01/KKS-05/UPS-001	
SERIAL NO : JS08T2L	SERIAL NO : 1526803063AF413	
MODEL : LENOVO X3500 M5	MODEL : EmersonGXTB-1000RT230	

If Alert is tick (X), fill in & elaborate details in FINDINGS / COMMENTS / DETECTION session.

PART I - HARDWARE, CABLING & CONNECTIVITY (INSPECTION & PHYSICAL CHECK)

Standard Task:	SERVER 1		SERVER 2		UPS	
	Function	Not Function & Resolved	Function	Not Function & Unresolved	Function	Not Function & Unresolved
Cleaning	<input checked="" type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>	
Verify Condition :						
Connection via LAN / USB	<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>	
Fan/ Power Supply	<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>	
Light Path Diagnosis panel	<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>	
UTP Cable connect to switch	<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>	
HA Gateway to Core Switch	<input type="checkbox"/>		<input type="checkbox"/>		<input type="checkbox"/>	
<i>(For Hospital only)</i>						

By pass Test	<input type="checkbox"/>	Alert	<input type="checkbox"/>	Not Function & Unresolved	
Webboard	<input type="checkbox"/>	Alert	<input type="checkbox"/>	Not Function & Unresolved	
Multilink	<input type="checkbox"/>	Alert	<input type="checkbox"/>	Not Function & Unresolved	

Backup Tape	Available	Type	Qty
	5	RDX	

PART II - SETUP & HEALTH CHECK

Function	Not Function & Resolved	Not Function & Unresolved
Run DSA & collect DSA Log		
SERVER 1	<input type="checkbox"/>	Alert
SERVER 2	<input type="checkbox"/>	Alert

Software/Firmware Version :	SERVER 1	SERVER 2
ClamAV	Current : 0.99.1	Updated
Zabbix	2.2.6	
Bacula	5.2.12	
Glassfish	3.1.2.2	
UEFI	1.11	
IMM2	1.75	

SERVER 1	SERVER 2
Current : Updated	Current : Updated
SUSE 11	
3-0.76-0.11	
4.3.4	

SERVER 1	SERVER 2
Current : Updated	Current : Updated
Screenshot / attachment	Screenshot / attachment

System Software
SUSE LINUX (OS)
Kernel
GCC

Housekeeping:
 Clear /tmp folder > 2 weeks
 Archive Log Files

N/A

**PREVENTIVE MAINTENANCE CHECKLIST
(SERVER / UPS / STORAGE)**

(FORM 2)

PAT ID : _____
 SERIAL NO: _____
 TYPE : STORAGE

N/A

Function	Not Function & Resolved		Not Function & Unresolved Alert	Ticket No
	Alert	Ticket No		
Perform Health Check:				
Storage Manager Error Alerts	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Storage events logs	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
HDD indicator lights	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Enable read/write cache	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Enable Media Scan both at controller level and LUN level	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Software/Firmware Version :

Controller	Updated
ESM	
Driver	

PART IV - FINDINGS / ALERTS / DETECTION

PART V - ACKNOWLEDGEMENT

PHARMAIAGA LOGISTICS SDN BHD (PLSB)

Signature & Stamp Customer Name Date	Signature & Stamp Engineer Name Date
 NOOR SHAहरुUL BIN MUSA RPh 6411 Pegawai Farmasi Klinik Kesihatan Kuala Perlis	 Ahmad Redzuan Bin Aziz Field Technical Engineer
28/2/19	28/02/2016

**PREVENTIVE MAINTENANCE CHECKLISTS
(CORE SWITCH / DISTRIBUTION SWITCH)**

(FORM 3)

FACILITY : KLINIK KESIHATAN KUALA PERLIS DEPT : IT STATE : PERLIS ZONE : 1

PAT ID :	CPS/PLS-01/KKS-05/DSW-001	LOCATION :	Location :	Outpatient Pharmacy
SERIAL NO:	DN120160414	Loaned <input type="checkbox"/>	Unlocated <input type="checkbox"/>	Relocated <input type="checkbox"/>

PART I - STANDARD TASK					PART II - STANDARD HEALTH CHECK					PART III - STANDARD SETUP				
Function	Not Function & Resolved	Alert	Not Function & Unresolved	Ticket No	Function	Not Function & Resolved	Alert	Not Function & Unresolved	Ticket No	Function	Not Function & Resolved	Alert	Not Function & Unresolved	Ticket No
Cleaning <input checked="" type="checkbox"/>	Verify Condition : Cable & Connection	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		SSH Access & Connection	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Execute Putty to start	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Power Supply	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Specify the IP to connect	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Test Login	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Check LED Status	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Test & Verify Management IP via console access	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	System LED	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>										
Master LED	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>											
STACK ID	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>											

PAT ID : _____ LOCATION : _____ Location : _____

SERIAL NO: _____ Loaned Unlocated Relocated

PART I - STANDARD TASK					PART II - STANDARD HEALTH CHECK					PART III - STANDARD SETUP				
Function	Not Function & Resolved	Alert	Not Function & Unresolved	Ticket No	Function	Not Function & Resolved	Alert	Not Function & Unresolved	Ticket No	Function	Not Function & Resolved	Alert	Not Function & Unresolved	Ticket No
Cleaning <input type="checkbox"/>	Verify Condition : Cable & Connection	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		SSH Access & Connection	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Execute Putty to start	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Power Supply	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Specify the IP to connect	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Test Login	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Check LED Status	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Test & Verify Management IP via console access	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	System LED	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>										
Master LED	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>											
STACK ID	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>											

PAT ID : _____ LOCATION : _____ Location : _____

SERIAL NO: _____ Loaned Unlocated Relocated

PART IV - FINDINGS / COMMENTS / DETECTION
PART V - ACKNOWLEDGEMENT

CUSTOMER Signature & Stamp Customer Name Date	PHARMANIAGA LOGISTICS SDN BHD (PLSB) Signature & Stamp Engineer Name Date
---	--



Appendix C

Hardware Inventory List



PHIS & CPS PROJECT

pharmaniaga

Hardware Inventory List

List of equipment available in the facility as per below table:

1. PCD

No	CPU	MONITOR	KEYBOARD	MOUSE	ASSET TAGGING	DEPARTMENT	LOCATION
1	PC05FLER	V901TZT4	0002536	HS515HA1K2Z	CPS/PLS-01/KKS-05/PCD-001	Outpatient Pharmacy	Screening, Verification, Labelling & Dispensing
2	PC04MX7U	V9022NMW	01282950	HS509HAOK38	CPS/PLS-01/KKS-05/PCD-002	Outpatient Pharmacy	Screening, Verification, Labelling & Dispensing 1
3	PC05FL1R	V902D1F8	0005120	44F5738	CPS/PLS-01/KKS-05/PCD-003	Outpatient Pharmacy	Medication Counselling & MTAC

2. THERMAL PRINTER

No.	THERMAL S/N	ASSET TAGGING	DEPARTMENT	LOCATION
1	2305Y001491	CPS/PLS-01/KKS-05/THP-001	Outpatient Pharmacy	Screening, Verification, Labelling & Dispensing 1

3. LASERJET PRINTER

No.	LASERJET S/N	ASSET TAGGING	DEPARTMENT	LOCATION
1	45146PLM2X2BR	CPS/PLS-01/KKS-05/LJP-001	Outpatient Pharmacy	Screening, Verification, Labelling & Dispensing 1

4.NETWORK SWITCH

No.	SWITCH S/N	ASSET TAGGING	DEPARTMENT	LOCATION
1	DN120160414	CPS/PLS-01/KKS-05/DSW-001	IT	Outpatient Pharmacy

5.SERVER

No.	SERVER S/N	ASSET TAGGING	DEPARTMENT	LOCATION
1	J308T2L	CPS/PLS-01/KKS-05/SVR-001	IT	Outpatient Pharmacy

6.UPS

No.	UPS S/N	ASSET TAGGING	DEPARTMENT	LOCATION
1	1526803063AF413	CPS/PLS-01/KKS-05/UPS-001	IT	Outpatient Pharmacy