












**Proposed Solution**

No.	Ticket #	Description	Requests	Proposed Solution
1.	I-PHIS021514717S	MTAC reporting - Request adherence score non mandatory	To remove mandatory for Adherence Score field.	To remove mandatory for Adherence Score field
2.	I-PHIS021577017S	Here is the issue: We need to cancel the wrong scheduled appointment date 19/7/17 (tarikh yang sudah lepas) by today. However, it ask us to make the patient visit (in actual, patient does not show up) and system ask us to change the appointment date (which is not right either as patient original appointment to be cancelled is 19/7/17 and not today).	To link MTAC record with previous visit if patient already have visit registered previously	To link MTAC record with previous visit if patient already have visit registered previously
3.	I-PHIS017801417S	Request by user for MTAC order refer below :  Type Of MTAC = Warfarin For button "Search History" to show Target INR Range and Bleeding Risk Factor of last visit at MTAC History only for Warfarin only.	Request by user for MTAC order refer below :  Type Of MTAC = Warfarin  For button "Search History" to show Target INR Range and Bleeding Risk Factor of last visit at MTAC History only for Warfarin only.	1) Search History button - to display Target INR Range and Bleeding Risk Factor of last visit at MTAC History only for Warfarin 2) Warfarin dose recommendation to be add in the MTAC Warfarin specific details sections

Function Flow	MTAC Screen in current Version 1.7	MTAC Screen in new Version 1.8
<p><b>1. MTAC Reporting – Diabetes Mellitus</b></p>	<p><b>Pharmacy Transaction → Clinical → MTAC Reporting → Diabetes Mellitus</b></p> <p>Currently there are mandatory to fill in Adherence Score for MTAC – Diabetes Mellitus</p> <p>1) Open Reporting section, click on Save without fill in the Adherence Score</p>  <p>2) Click on Save</p> 	<p><b>Pharmacy Transaction → Clinical → MTAC Reporting → Diabetes Mellitus</b></p> <p>To remove mandatory for Adherence Score field.</p> <p>1) Open Reporting section, click on Save without fill in the Adherence Score</p>  <p>2) Click on Save</p> 

Function Flow	MTAC Screen in current Version 1.7	MTAC Screen in new Version 1.8
<p><b>1. MTAC Reporting – Visit</b></p>	<p><b>Pharmacy Transaction → Clinical → MTAC Reporting</b></p> <p>Currently User need to register new visit if patient don't have active visit to proceed with MTAC Reporting</p> <ol style="list-style-type: none"> <li>1) Open MTAC Reporting listing screen</li> <li>2) Select patient with visit un-active</li> <li>3) Pop-up message will appear – Please record patient visit first</li> </ol> 	<p><b>Pharmacy Transaction → Clinical → MTAC Reporting</b></p> <p>To link MTAC record with previous visit if patient already have visit registered previously</p> <ol style="list-style-type: none"> <li>1) Open MTAC Reporting listing screen</li> <li>2) Select patient with visit un-active</li> <li>3) System will refer to patient previous visit if any and will allow to proceed with MTAC Reporting screen</li> </ol>   <p><i>Note: If patient is first time and dont have any visit at all, user may need to register the visit first.</i></p>

Function Flow	MTAC Screen in current Version 1.7	MTAC Screen in new Version 1.8
<p><b>1. MTAC – Warfarin</b></p>	<p><b>Pharmacy Transaction → Clinical → MTAC Order</b></p> <p>1) Search History button - to display Target INR Range and Bleeding Risk Factor of last visit at MTAC History only for Warfarin</p> <p>Step:</p> <ol style="list-style-type: none"> <li>Open MTAC Order Screen</li> <li>Select Type of MTAC = Warfarin</li> </ol> <p><i>Note: Patient should have MTAC Warfarin History in system</i></p>  <p>c. Click on 'Search History' Button</p> <ul style="list-style-type: none"> <li>Currently No 'Target INR Range and Bleeding Risk Factor' been displayed</li> </ul> 	<p><b>Pharmacy Transaction → Clinical → MTAC Order</b></p> <p>1) Search History button - to display Target INR Range and Bleeding Risk Factor of last visit at MTAC History only for Warfarin</p> <p>Step:</p> <ol style="list-style-type: none"> <li>Open MTAC Order Screen</li> <li>Select Type of MTAC = Warfarin</li> </ol> <p><i>Note: Patient should have MTAC Warfarin History in system</i></p>  <p>c. Click on 'Search History' Button</p> <ul style="list-style-type: none"> <li>'Target INR Range and Bleeding Risk Factor' will be displayed</li> </ul> 

Function Flow	MTAC Screen in current Version 1.7	MTAC Screen in new Version 1.8
	<p><b>Pharmacy Transaction → Clinical → MTAC Reporting → Specific Details</b></p> <p>2) Warfarin dose recommendation to be add in the MTAC Warfarin specific details sections</p> <p>Currently No MTAC Warfarin specific details sections in Specific Details MTAC Warfarin Screen</p> 	<p><b>Pharmacy Transaction → Clinical → MTAC Reporting → Specific Details</b></p> <p>2) Warfarin dose recommendation has been added in the MTAC Warfarin specific details sections</p> 