














### Summary of Ticket Number and Requests for Radiopharmaceuticals - General

No.	Ticket #	Description	Requests	Solution
1.	I-PHIS013456716S	<b><u>Drug Master - Default dose/UOM for all RP screen</u></b> Default dose/uom - to apply in all screen, not only in order	To allow set default dose for all screens	1. Default dose/UOM for all RP screen
2.	I-PHIS013453116S	<b><u>Therapy Linked Radiopharma - Default setting</u></b> To set default in Facility	To set default in Facility	2. Radiopharma Therapy (Maintenance) - To remove default function in IWP and add function for facility user able to set default radiopharmaceuticals.

Function Flow	Drug Master – Default UOM in current Version 1.7	Drug Master – Default UOM in new Version 1.8
<p>1. Drug Master – Default UOM</p>	<ul style="list-style-type: none"> <li>If user set default UOM in drug master:            </li> <li>Impact to only Radiopharmaceuticals Appointment, Record Prescription and Medication Order screens below:             </li> </ul>	<ul style="list-style-type: none"> <li>New batch allocation flow:            </li> <li>Added more screens to refer same default Dose/UOM value:           <ol style="list-style-type: none"> <li>Calculate Activity button in Preparation screen                </li> <li>Generator Elution screen                </li> </ol> </li> </ul>

Function Flow	Drug Master – Default UOM in current Version 1.7	Drug Master – Default UOM in new Version 1.8																																																	
		<p>3. Kit Preparation screen</p>  <p>4. Radioisotope Receipt Form report</p>  <table border="1"> <thead> <tr> <th rowspan="2">Radiopharmaceutical/Product</th> <th rowspan="2">Brand Name</th> <th rowspan="2">Batch No.</th> <th colspan="3">Received</th> <th rowspan="2">Supplier</th> <th rowspan="2">Expiry Date/Time</th> <th colspan="2">Calibration</th> <th colspan="2">Measured</th> <th rowspan="2">Status</th> </tr> <tr> <th>Qty</th> <th>UOM</th> <th>Date</th> <th>Date</th> <th>Activity</th> <th>UOM</th> <th>Activity</th> <th>UOM</th> </tr> </thead> <tbody> <tr> <td>Sodium Perchnetate (Tc-99m)</td> <td>MOH</td> <td>9909</td> <td>10000</td> <td>vial</td> <td>02/04/2018</td> <td></td> <td>01/05/2018 03:05 AM</td> <td>05/04/2018 03:05 PM</td> <td>67.57</td> <td>mCi</td> <td>63.93</td> <td>mCi</td> <td></td> </tr> <tr> <td>Sodium Perchnetate (Tc-99m)</td> <td>MOH</td> <td>G3454</td> <td>1000</td> <td>vial</td> <td>09/03/2018</td> <td></td> <td>30/04/2018 11:04 AM</td> <td>09/03/2018 11:03 AM</td> <td>27.03</td> <td>mCi</td> <td>27.03</td> <td>mCi</td> <td></td> </tr> </tbody> </table>	Radiopharmaceutical/Product	Brand Name	Batch No.	Received			Supplier	Expiry Date/Time	Calibration		Measured		Status	Qty	UOM	Date	Date	Activity	UOM	Activity	UOM	Sodium Perchnetate (Tc-99m)	MOH	9909	10000	vial	02/04/2018		01/05/2018 03:05 AM	05/04/2018 03:05 PM	67.57	mCi	63.93	mCi		Sodium Perchnetate (Tc-99m)	MOH	G3454	1000	vial	09/03/2018		30/04/2018 11:04 AM	09/03/2018 11:03 AM	27.03	mCi	27.03	mCi	
Radiopharmaceutical/Product	Brand Name	Batch No.				Received					Supplier	Expiry Date/Time	Calibration			Measured		Status																																	
			Qty	UOM	Date	Date	Activity	UOM	Activity	UOM																																									
Sodium Perchnetate (Tc-99m)	MOH	9909	10000	vial	02/04/2018		01/05/2018 03:05 AM	05/04/2018 03:05 PM	67.57	mCi	63.93	mCi																																							
Sodium Perchnetate (Tc-99m)	MOH	G3454	1000	vial	09/03/2018		30/04/2018 11:04 AM	09/03/2018 11:03 AM	27.03	mCi	27.03	mCi																																							

Function Flow	Drug Master – Default UOM in current Version 1.7	Drug Master – Default UOM in new Version 1.8
<p>2. Radiopharma Therapy</p>	<ul style="list-style-type: none"> <li>Current screen – IWP: HQ user able to set default</li> </ul>  <ul style="list-style-type: none"> <li>Current screen - Facility: Facility user unable to set the default column, value is set in HQ only</li> </ul> 	<ul style="list-style-type: none"> <li>Remove default function in HQ; and to allow facility to set default Radiopharmaceuticals:</li> </ul>  <ul style="list-style-type: none"> <li>Allow facility to set default Radiopharmaceuticals:</li> </ul> 